

MODEL CR-79 CAGE WASH RACK



MODEL CR-79 CAGE WASH RACK

LYNX Universal Wash Racks are versatile, heavy-duty, stainless steel racks designed to securely hold rodent cages, filter tops, and debris pans for thorough, efficient cleaning in a walk-in type Cage and Rack Washer.

STANDARD FEATURES:

- Adjustable Retainers
- Adjustable Shelf Positioning
- Third Wheel Provides Zero Degree Turn
- Lightweight for Easy Maneuvering
- Four (4) Pair of Cage Wash Racks
- Ultra High Capacity (consult factory)

CONSTRUCTION:

- Complete 304 Stainless Steel Construction
- Six (6) Stainless Steel Casters
- Polyolefin Wheels with Stainless Steel Axles and Delrin Bearings

OPTIONAL FEATURES:

- Additional Cage Wash Rack Shelves
- Additional Pan Wash Rack Shelves
- Additional General Purpose Rack Shelves

Size:

Overall Unit:

34³/₈" wide x 82¹/₂" high x 84" long


Additional Models: CR-69, CR-39

Note: Wash Racks can be designed to any size to meet new or existing applications.

ACCESSORIES


CR-79 CAGE WASH RACK

CR-79	83H x 35W x 84L
CR-69	83H x 35W x 74L
CR-39	83H x 35W x 44L



Adjustable retainer bar to accommodate a variety of cage, filtertop, debris pan, and accessory sizes.

Adjustable Shelving Position



CR-79 Option Shown Here



For further information, please contact:

CAGE WASH RACKS

