


The LYNX Advantage...

Bedding Handling Systems

Why you should choose LYNX for your Waste Bedding Disposal and Clean Bedding Delivery Systems...

- *Bedding Solutions Customized for your Facility!*
- *Advanced Ergonomics and Safety*
- *High Throughput and Efficiency*
- *Extremely Cost Effective*
- *Will Handle Most Bedding* Types*
- *Engineered for Easy Maintenance*
- *Utilizes LYNX Smart Purge™ System*
- *Highly Advanced Vacuum Technology*

*Equipment designed for free-flowing and non-bridging bedding




LYNX is Committed to:

- Use of Non-Proprietary Parts
- Exceptional Customer Service
- State-of-the-art Microprocessor Program
- Rugged, Reliable Construction



Typical Waste Bedding Disposal Components...

- 304 S/S Waste Bedding Dump Station (stand alone or integral)
- Engineered 304 S/S Pipe Network to prevent clogging
- Waste Collection Module
 - 304 S/S Waste Bedding Intercept Hopper
 - Sealed, Vacuum Rated Waste Container
- 99.9% Efficient Primary Filter Separator
- Allergen Dust Collector (optional)
- Odor Control Filter
- Central Vacuum Power Unit
- System Controls Package



Waste Bedding
Dump Station



Vacuum Rated
Waste Container




Allergen Dust
Collector (optional)



Primary Filter
Separator

Typical Clean Bedding Delivery Components...

- Bulk Bag Discharge Station
- Engineered 304 S/S Pipe Network to prevent Clogging
- 304 S/S Filtered Material Receiver
- HEPA Intake Filter for Ultra Fine Dust Capture
- Central Vacuum Power Unit
- System Controls Package
- Bedding Dispenser (or Retrofit to Existing Dispenser)



Material Receiver



Central Vacuum Pump
Unit



Bulk Bag Discharge
Station

